

DONJOY®

Lumbostrap 35
Dorsostrap 55

SUPPORT AND COMFORT REFINED



enovis™



DONJOY®

Support and comfort refined

While effective support should be standard, comfort is sometimes an afterthought. With the new DonJoy Lumbostrap 35 and DonJoy Dorsostrip 55 back braces, you can have both. High-quality materials and semi-rigid construction help deliver new levels of comfort for patients requiring support and stability for the low back.





DonJoy Lumbostrap 35

Semi-rigid lumbar brace with a height of 35 cm, for support and stabilization of the lumbar spine from S1-T10 and relief of low back pain.

INDICATIONS

- Results of vertebral collapse (most common indication), that requires a very long usage (up to 12 months)
- Outcomes of vertebral fractures of the lower dorso-lumbar tract
- Post-operation period
- Lumbalgia
- Sciatic lumbalgia and lumbar cruralgia
- Spondylarthrosis and lumbar discopathies
- Slight traumas to the lumbosacral column and paravertebral lumbar contractures
- Treatment of fracture of transverse processes
- Degenerative lumbar scoliosis
- Spondylosis with or without vertebral listhesis

FEATURES AND BENEFITS

- Four semi-rigid pre-shaped **dorsal stays** help provide adjustable enhanced back support
- Two of the stays can be replaced with a **more rigid** pair for additional stability
- **Adjustable straps** provide individual compression
- **Anatomically-contoured** for improved comfort
- **Lightweight, breathable, odor-resistant materials** aid hygiene and compliance



ORDERING INFORMATION

BRACE HEIGHT: 35cm

PART NUMBER	WAIST CIRCUMFERENCE	SIZE
82-0055-1	56 - 68 cm (22 - 26 1/2")	XS
82-0055-2	68 - 82 cm (26 1/2 - 32 1/4")	S
82-0055-3	82 - 98 cm (32 1/4 - 38 1/2")	M
82-0055-4	98 - 116 cm (38 1/2 - 45 1/2")	L
82-0055-5	116 - 136 cm (45 1/2 - 53 1/2")	XL
82-0055-6	136 - 156 cm (53 1/2 - 61 1/4")	XXL





DonJoy Dorsostrap 55

Semi-rigid dorso-lumbar brace with shoulder straps and a height of 55 cm. For support and stabilization of the dorsal-lumbar spine from S1-T4 and relief of low back pain.

INDICATIONS

- Outcomes of vertebral fractures of the dorso-lumbar tract
- Dorsal vertebral collapses on an osteoporotic or metastatic basis (most common indication), that requires a very long usage (up to 12 months)
- Post-operation period
- Osteoporosis
- Osteomalacia
- Dorsal kyphosis
- Osteolysis
- Kyphosis with chronic back pain
- Muscular insufficiency

FEATURES AND BENEFITS

- **Reverse shoulder** straps to help improve kyphotic posture
- Four semi-rigid pre-shaped **dorsal stays** help provide adjustable enhanced back support
- Two of the stays can be replaced with a more **rigid pair** for additional stability
- **Adjustable straps** provide individual compression
- **Anatomically-contoured** for improved comfort
- **Lightweight, breathable, odor-resistant materials** aid hygiene and compliance



ORDERING INFORMATION

BRACE HEIGHT: 55cm

PART NUMBER	WAIST CIRCUMFERENCE	SIZE
82-0056-1	56 - 68 cm (22 - 26 1/2")	XS
82-0056-2	68 - 82 cm (26 1/2 - 32 1/4")	S
82-0056-3	82 - 98 cm (32 1/4 - 38 1/2")	M
82-0056-4	98 - 116 cm (38 1/2 - 45 1/2")	L
82-0056-5	116 - 136 cm (45 1/2 - 53 1/2")	XL
82-0056-6	136 - 156 cm (53 1/2 - 61 1/4")	XXL

DONJOY®

Asia Pacific

1905, Tower 2, Grand Central Plaza
Shatin
HONG KONG
T: 852 31051415
E: marketing.ap@enovis.com

Australia – New Zealand

PO Box 6057
Frenchs Forest DC
NSW 2086
AUSTRALIA
T: 1300 66 77 30
E: customerservice.au@enovis.com

Benelux

Kleinhoefstraat 5, bus 39
Geel 2440
BELGIUM
T: +32 (0)14 248090
E: orders.bx@enovis.com

Canada

6485 Kennedy Rd
Mississauga
Ontario, L5T2W4
CANADA
T: 1-866-866-5031
E: canada.orders@enovis.com

China

A312 SOHO ZhongShan Plaza
1055 W. ZhongShan Rd
Shanghai, 200051
CHINA
T: (8621) 6031 9989
E: shujun.tang@enovis.com

France

3 rue de Bethar
Centre Européen de Frêt
64990 Mouguerre
FRANCE
T: +33 (0)5 59 52 86 90
E: sce.cial@enovis.com

Germany

Bötzingen Straße 90
79111 Freiburg
GERMANY
T: +49 (0) 180 1 676 333
E: kundenservice@enovis.com

Italy

Via Leonardo da Vinci 97
Trezzano sul Naviglio
Milan 20090
ITALY
T: +39 02 835 98001
E: vendite@enovis.com

Nordic

Murmansgatan 126
212 25 Malmö
SWEDEN
T: +46 40 39 40 00
E: info.nordic@enovis.com

South Africa

Unit 6, 2 on London, 2 London Circle
BrackenGate Business Park
Brackenfell, 7560
SOUTH AFRICA
T: +27 (0) 21 276 2968
E: info.southafrica@enovis.com

Spain

Av.Cornella 144. 1-4
08950 Esplugues de Llobregat
Barcelona
SPAIN
T: +34 93 480 32 02
E: ventas@enovis.com

UK & Ireland

1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UK
T: +44 (0)1483 459 659
E: ukorders@enovis.com

USA

2900 Lake Vista Drive
Lewisville
TX 75067
U.S.A.
T: 760 727 1280
E: customercare@enovis.com

enovis™

Creating Better Together™

enovis-medtech.eu

1a Guildford Business Park, Guildford, Surrey, GU2 8XG, UNITED KINGDOM

Individual results may vary. Neither Enovis, DJO, LLC or any of their subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

MKT-000-00090550-EN - Rev A

Copyright © 2022 DJO, LLC, a subsidiary of Enovis Corporation